

MHP-100 A/S Manual

☒ mat

☒ regulator

AS-20130401

Nuga Medical Customer Support Team


	division	standard	note
Indication	Total product length	580 (W) * 760 (D) * 15 (H)	
	Product weight	3.3kg / 1SET	
	Total product power consumption	MAX 150W ± 15%.	
	electrical rating	AC220V / 50Hz/ 60Hz	
	Manufacturing license number	No. 1478	
	Item License Number	10-769	
Specifications	thermic rays	220V/165W/320Ω±10%	
	temperature Sensor	10kΩ	
	Protector	100°C ± 5°C× 3	

MHP-100 Error Code and Action

mark	Abnormal contents	cause	Confirmation details	action
	Temperature sensor disconnection	Error due to non-detection of the mat's temperature sensor (temperature sensor system disconnection)	1. Controller-Mat: Check connection status 2. Regulator cable: Red-brown single wire (continuity test) 3. Disconnection of the temperature sensor system in the mat 4. Mat connector pin (5 pins): Check soldering area	1. Regulator-Mat: Reconnection 2. Connect broken wire or replace cable ASSY 3. Connect lead wire or replace temperature sensor 4. Check soldering area
	Temperature sensor short circuit	Temperature sensor cable CIRCUIT-SHORT	Temperature sensor defective	Temperature sensor replacement (requires sewing)
	Heating element overheating	When the temperature sensor detects a temperature above 86°C	1. Actual mat temperature overheating 2. In case of temperature detection error due to temperature sensor error	1. Check the cause of overheating and take action 2. Abnormal after measuring temperature sensor resistance value case replacement
 On the regulator Display temperature is normal > mat heating is poor	Poor fever	The regulator is normal Heat on the mat If it doesn't work	1. Hot wire disconnection 2. Bimetal lead wire disconnection 3. Bimetal defect (low occurrence rate)	1. After removing the mat, visually check the hot wire breakage area >> Replace the hot wire ASSY 2. Separate the mat and check for disconnection >> Solder the lead wire (connection) 3. Replace the bimetal


How to check when the mat has poor heating

Temperature sensor system regeneration repair
method Bimetal system regeneration repair method



Malfunction symptoms	Checklist	Action taken
<p>Even after pressing the power button</p> <p>The controller does not power on.</p> <p>The power lamp (green) does not light up.</p> 	<ul style="list-style-type: none"> - Check the power plug and outlet connection status. 	<ul style="list-style-type: none"> - Plug the power cord into the outlet correctly.
	<ul style="list-style-type: none"> - Check fuse disconnection 	<ul style="list-style-type: none"> - Replace the fuse. See page 6 * If the fuse continues to disconnect even after replacement, replace it with another mat and test. If the same symptom occurs, replace the PCB.
	<ul style="list-style-type: none"> - Check AC110/220V measurement on terminals 1 and 2 See page 6 	<p>Power cable disconnection</p> <p>→ Replace cable A'SSY.</p>
	<ul style="list-style-type: none"> - Check DC5V measurement across ZD1 	<ul style="list-style-type: none"> - If 5V is measured, replace the power LED. (54, RED) - If measurement is not performed, replace the PCB.
<p>The mat does not work.</p> <p>(Heater LED does not light up.)</p>	<ul style="list-style-type: none"> - Check the connection status of the mat <p>(If the sensor is disconnected, the heater does not operate.)</p>	<ul style="list-style-type: none"> - Test by replacing the cable and see if the heater works. If not, replace a mat that works normally and test it. I'll try it. If it does not operate properly after replacing the mat, replace the PCB.
<p>Temperature control is not possible.</p>	<ul style="list-style-type: none"> - Check variable resistance status <p>(Check for damage in one direction check if it keeps running)</p>	<ul style="list-style-type: none"> - Replace the variable resistor. - Replace the PCB.

MHP100 controller PCB main parts



1-2: 110/220V power input
 3-4: Heated wire power supply
 5-6: Temperature sensor

- ① Time setting Dip Switch: 6 to 12 hour settings.
- ② Power Tact Switch: ITS1105(H 8mm)
- ③ AC 220V output control TRIAC: BTA 06 6000B
- ④ Heating Pad temperature control V-Resistor: 10kΩ, 160, 210
- ⑤ Power FUSE: 250V 1.5A
- ⑥ Overcurrent protection Varistor: 10D471 (100 470V)




Currently, the shape of the time control switch has been changed to 8H (left) and 12H (right).

① Time adjustment according to Dip Switch settings


S/W 1	S/W 2	Hour	
ON	ON	6Hr	
ON	OFF	8Hr	
OFF	ON	9Hr	
OFF	OFF	12Hr	


※ Set to 8 hours when the product is shipped.





number	Part Name	function	Malfunction symptoms	Diagnosis method	Action method
2	TACT S/W	ON/OFF	ON/OFF does not work well.	You have to press it hard to turn it on/off.	switch replacement
3	Triac	Output control according to temperature setting	Temperature control is not possible. (Example) > It is set to 40 degrees, but the temperature continues to rise. > The mat heats up even when the switch is turned off.	Even when the switch is OFF, the voltage is checked to be around 100V between terminals 1 to 4 of the connector.	Triac replacement
4	VR	Setting temperature control	The correct temperature setting is not possible.	Lightly tap the area around the temperature control knob (yellow dotted line in the picture above) to change the temperature.	VR cleaning (using high pressure air) or replacement
5	Fuse	Overcurrent blocking (2A)	The controller does not turn on.	Check for blown fuse in the controller PCB.	Fuse replacement

1. Heating Wire : AC 220V
150W \pm 15%.
297 Ω \pm 10%.
2. Protector : 100 $^{\circ}$ C \pm 5%
3. Temperature Sensor : 10k Ω







Temperature sensor resistance: around $10K\Omega$ (based on $25^{\circ}C$)
Heating resistance: around 2972



MHP100 Mat-Regulator Schematic




When the black or green cable of the controller cable is disconnected
→ the same phenomenon as symptom number 2 on page 4 appears.


When the red or brown cable of the controller cable is disconnected →
the same phenomenon as symptom number 1 on page 4 appears.





MHP100 regulator terminal layout




1	thermic rays
2	temperature Senger
3	temperature Senger
4	thermic rays
5	ground connection




MHP100 heated wire & bimetal lead wire wiring detailed diagram




MHP100 failure location







Symptom	<p>Before the disconnection part (mat 1/2) heated.</p> <p>Heating resistance value: Measured at approximately 150 yaw.</p>	<p>Before the disconnected part (3/4 of the mat) was heated: Heating wire resistance value: Measured at approximately 200 yaw.</p>	<p>Poor heat generation throughout the mat.</p>
cause	<p>1. Hot wire disconnection in part (1)</p> <p>► As the core wire and coil stick together (short circuit), In case of disconnection.</p>	<p>1. Hot wire disconnection in part (2)</p> <p>► When the core wire and coil are attached (short-circuited) and disconnected.</p>	<p>1. Disconnection of heating wire</p> <p>► When the coil and core wire are not attached</p> <p>2. Lead wire disconnection</p> <p>► Bimetal and temperature sensor lead wire</p> <p>3. Regulator failure</p>



Precautions when using MHP100 product (1)



Area where disconnection occurs

Make sure the mat is not folded at an acute angle as shown. (Cause of central hot wire breakage)

It is best to unfold and store the product after use!!!!



If you have to fold it...



Let the turbanium touch each other...



of the

Fold the folded part so that it forms a round shape.

Precautions when using MHP100 product (2)



How to replace MHP-100 heated ASSY



mat



◀ To replace the heating wire, separate it into the bottom plate, heating wire ASSY, and turbanium fabric (3 PCS).


Place the heated ASSY to be replaced on the bottom plate. ▶





◀ It is treated with hot melt to attach the bottom plate and hot wire ASSY.

◀ Apply hot melt to the inside of the seam along the border of the product.



Process inside the seam line
along the border of the product.





◀ Cut the scrap portion of the heated ASSY to fit the fabric.


The size of the heated ASSY is slightly larger than that of the existing product, so it must be cut.

→ At this time, be careful not to damage the heating and lead wires!!

Be careful not to damage the heating and lead wires!!

Apply hot melt to the inside of the seam along the border of the product. ▼

Process inside the seam line along the border of the product.



1. When the upper (center) part of the mat's heating wire is disconnected.




1. Pull the left/right heating wire toward the center to repair the disconnected part of the heating wire.




For the left part, pull the heating wire toward the center in the same way.




1. Pull the left/right heating wire toward the center to repair the disconnected part of the heating wire.



1. Hot wire - How to connect hot wire





hot wire to connect



Carefully peel off the red silicone outer covering and separate the coil heating wire and core wire.


1. Hot wire - hot wire connection details




Before connecting the hot wire, insert the shrink tube in advance.

1. Detailed diagram of hot wire (core wire-core wire) connection


To connect the core wires, use weak soldering or twist them with a twist, and then pull the shrink tube on both sides to complete the double closure.



1. Detailed diagram of heating wire (COIL-COIL) connection



1. Finishing the hot wire connection



1. Glass fiber finishing

fiber glass
Finish with
(silver tape)



1. After finishing, measure the resistance value of the hot wire

