

NUGA PRODUCT A/S GUIDE

SMALL PRODUCT

NM-80 / (TP80)




NM-80




HEATING WIRE ASSY / COIL ASSY


WIRE CONNECTION: SERIES TYPE



PROBLEM	PARTS	SOLUTION
NO HEAT	WIRE DISCONNECTION/ BIMETAL WIRE/SENSOR/ REMOTE PCB/HEATING WIRE ASSY	REPLACE/REPAIR
RANDOM TEMPERATURE DISPLAY	BROKEN SENSOR / DISCONNECT SENSOR WIRE	REPLACE SENSOR/ REPLACE POTENTIOMETER/VOLUME METER
REMOTE SUDDENLY NO POWER AFTER INSERTING TO HEATING MAT / FUSE ALWAYS EXPLODE	TERMINAL PIN	REPLACE TERMINAL PIN REPLACE FUSE




OPENING THE BODY




NM-80



REMOTE CONTROL



OLD MODEL



NEW MODEL

DISPLAY ERROR



- TEMPERATURE SENSOR (DISCONNECT / BROKEN)



- TEMPERATURE SENSOR / REMOTE CONTROL PCB

REMOTE PCB PARTS

A – SWITCH

B – FUSE

C – TRANSISTOR

D – RESISTOR

E – CAPACITOR

F – POTENTIOMETER/VOLUME

PROBLEM	PARTS	SOLUTION
NO HEAT	REMOTE PCB/MAT	REPLACE PCB / SENSOR/ BIMETAL
NO POWER	REMOTE CABLE/FUSE/ SWITCH / CAPACITOR /RESISTOR	REPLACE REMOTE CABLE/FUSE/ SWITCH / CAPACITOR /RESISTOR / PROJECTOR CABLE
OVERHEAT	TRANSISTOR	REPLACE TRANSISTOR
DISPLAY RANDOM	POTENTIOMETER/VOLUME	REPLACE POTENTIOMETER / VOLUME / CHECK CONTACT TO PCB

NM-90




NM-90



HEATING WIRE ASSY / COIL ASSY

WIRE CONNECTION: SERIES TYPE




USUALLY DISCONNECTION OCCUR





PROBLEM	PARTS	SOLUTION
NO HEAT	WIRE DISCONNECTION/ BIMETAL WIRE/SENSOR/ REMOTE PCB/HEATING WIRE ASSY	REPAIR WIRE DISCONNECTION. REPLACE BIMETAL WIRE/SENSOR/ REMOTE PCB/HEATING WIRE ASSY
RANDOM TEMPERATURE DISPLAY	BROKEN SENSOR / DISCONNECT SENSOR WIRE	REPLACE / RECONNECTION
REMOTE SUDDENLY NO POWER AFTER INSERTING TO HEATING MAT / FUSE ALWAYS EXPLODE	TERMINAL PIN	REPLACE TERMINAL PIN REPLACE FUSE


OPENING THE BODY



BIMETAL




NM-90



OLD MODEL

REMOTE CONTROL



NEW MODEL

DISPLAY ERROR



- TEMPERATURE SENSOR (DISCONNECT / BROKEN)



- TEMPERATURE SENSOR / REMOTE CONTROL PCB

REMOTE PCB PARTS

A – SWITCH

B – FUSE

C – TRANSISTOR

D – RESISTOR

E – CAPACITOR

F – POTENTIOMETER/VOLUME

PROBLEM	PARTS	SOLUTION
NO HEAT	REMOTE PCB	REPLACE PCB / REMOTE ASSEMBLY
NO POWER	REMOTE CABLE/FUSE/ SWITCH / CAPACITOR /RESISTOR / PROJECTOR CABLE	REPLACE REMOTE CABLE/FUSE/ SWITCH / CAPACITOR /RESISTOR / PROJECTOR CABLE
OVERHEAT	TRANSISTOR	REPLACE TRANSISTOR
DISPLAY RANDOM	POTENTIOMETER/VOLUME	REPLACE POTENTIOMETER / VOLUME / CHECK CONTACT TO PCB


MHP-100



MHP-100



MHP100 열선, 온도센서, 바이메탈 배치도



1. Heating Wire : AC 220V
 $150W \pm 15\%$.
 $297\Omega \pm 10\%$.
2. Protector : $100^\circ\text{C} \pm 5\%$
3. Temperature Sensor : $10\text{k}\Omega$

PROBLEM	PARTS	SOLUTION
NO HEAT	BIMETAL WIRE	RECONNECTION
NO HEAT / RANDOM TEMPERATURE DISPLAY	BROKEN SENSOR / DISCONNECT SENSOR WIRE	REPLACE / RECONNECTION
NO HEAT	COIL WIRE	RECONNECT / REPLACE WHOLE COIL ASSEMBLY
REMOTE SUDDEN NO POWER AFTER INSERTING TO HEATING MAT / FUSE ALWAYS EXPLODE	TERMINAL PIN	REPLACE TERMINAL PIN REPLACE FUSE



MHP-100



MOST DISCONNECT PART


SYMBOL




MHP100 고장발생부위




고장발생부위




MHP100 조절기 단자 배치도



1	열선
2	온도센서
3	온도센서
4	열선
5	접지




MHP100 열선& 바이메탈 리드선 결선 상세도




MHP-100

MHP100 조절기 단자 배치도



1	열선
2	온도센서
3	온도센서
4	열선
5	접지




열선
온도센서
온도센서
열선
접지

MHP100 열선& 바이메탈 리드선 결선 상세도

Silicon coat



Teflon




Connection with core wire

심선과 연결

코일과 연결

Connection with coil

pvc



MHP-100



REMOTE CONTROL

REMOTE PCB PARTS

A – SWITCH


B – FUSE

C – TRANSISTOR

D – RESISTOR

E – CAPACITOR

F – POTENTIOMETER/VOLUME



DISPLAY ERROR




- TEMPERATURE SENSOR (DISCONNECT / BROKEN)



- TEMPERATURE SENSOR / REMOTE CONTROL PCB

PROBLEM	PARTS	SOLUTION
NO HEAT	REMOTE PCB	REPLACE PCB / REMOTE ASSEMBLY
NO POWER	REMOTE CABLE/FUSE/ SWITCH / CAPACITOR /RESISTOR /	REPLACE REMOTE CABLE /FUSE/ SWITCH / CAPACITOR /RESISTOR /
OVERHEAT	TRANSISTOR	REPLACE TRANSISTOR
DISPLAY RANDOM	POTENTIOMETER/VOLUME	REPLACE POTENTIOMETER / VOLUME / CHECK CONTACT TO PCB


T11






T11

OLD ALUMINUM COIL

PROBLEM	PARTS	SOLUTION
NO HEAT	REMOTE PCB / ALUMINUM COIL / SENSOR / BIMETAL	REPLACE REMOTE PCB / ALUMINUM COIL / SENSOR / BIMETAL
NO POWER	CHECK REMOTE CABLE FUSE/ SWITCH / CAPACITOR /RESISTOR	REPLACE REMOTE CABLE / FUSE/ SWITCH / CAPACITOR /RESISTOR
OVERHEAT	TRANSISTOR	REPLACE TRANSISTOR
DISPLAY RANDOM	POTENTIOMETER/VOLUME	REPLACE POTENTIOMETER / VOLUME / CHECK CONTACT TO PCB



BIMETAL




T9



T9

MAIN PCB



REMOTE PCB PARTS

A – FUSE

B – CAPACITOR (ORANGE)

C – TRANSISTOR

D – RESISTOR

E – CAPACITOR (BLACK)




MEMBRANE PCB



PROBLEM	PARTS	SOLUTION
NO HEAT	REMOTE PCB / ALUMINUM COIL/SENSOR/BIMETAL	REPLACE PCB / REMOTE PCB / ALUMINUM COIL/SENSOR/BIMETAL
NO POWER	POWER CORD / FUSE/ SWITCH / CAPACITOR /RESISTOR / PROJECTOR CABLE /MEMBRANE PCB	REPLACE POWER CORD /FUSE/ SWITCH / CAPACITOR /RESISTOR / PROJECTOR CABLE /MEMBRANE PCB
OVERHEAT	TRANSISTOR	REPLACE TRANSISTOR

T9

OPENING THE BODY




T5



T5

MAIN PCB



MAIN PCB PARTS



A – FUSE

B – TRANSISTOR

C – CAPACITOR (BLACK)

D – IC

E - RESISTOR CERAMIC



PROBLEM	PARTS	SOLUTION
NO HEAT	REMOTE PCB / ALUMINUM COIL/SENSOR/BIMETAL	REPLACE PCB / REMOTE PCB / ALUMINUM COIL/SENSOR/BIMETAL
NO POWER	CHECK POWER CORD/FUSE/ SWITCH / CAPACITOR /RESISTOR CERAMIC /IC	REPLACE POWER CORD/FUSE/ SWITCH / CAPACITOR /RESISTOR CERAMIC /IC
OVERHEAT	TRANSISTOR	REPLACE TRANSISTOR

NM-200



NM-200


ERROR CODE

Display	Descriptions	Display location
--	Disconnection of temperature sensor	The temperature sensor is not detected
E1	Short circuit of temperature sensor	The temperature sensor is under the short circuit
E2	Disconnection of heating unit	The heating unit is disconnected or under the disconnection by the bimetal operation
E3	Overheating of heating unit	The temperature in the heating unit is increased more than 86°C.

NM-200



REMOTE
PCB



REMOTE PCB PARTS

A – SWITCH

B – FUSE


C – TRANSISTOR

D – RESISTOR

E – CAPACITOR

F – POTENTIOMETER/VOLUME

PROBLEM	PARTS	SOLUTION
NO HEAT	TOURMANIUM HEATER / PROJECTOR CABLE / REMOTE PCB	REPLACE / DISCONNECT BROKEN WIRE OF PROJECTOR CABLE
NO POWER	REMOTE PCB	REPLACE FUSE/ SWITCH / CAPACITOR /RESISTOR / PROJECTOR CABLE
OVERHEAT	REMOTE PCB	REPLACE TRANSISTOR
DISPLAY RANDOM	POTENTIOMETER /VOLUME	REPLACE / CHECK CONTACT TO PCB



LAMP PCB



TOURMANIUM
HEATER



PROJECTOR
CABLE



TYPES OF
ERROR

MB2



MB2



MAIN PCB



BATTERY



ADAPTER/CHARGER

PROBLEM	PARTS	SOLUTION
NO POWER	SWITCH / BATTERY / MAIN PCB / CHARGER	REPLACE
NO FREQUENCY	MAIN PCB / PE CABLE	REPLACE

LF RUBBER SILICON (S)



LF RUBBER SILICON (B)



WASER



PE CABLE


NM300



NM300



MAIN PCB



PROBLEM	PARTS	NO.
NO POWER	FUSE	REPLACE 1
FUSE ALWAYS EXPLODE / ALL BULB EXPLODE EVEN NEW / NO POWER	ic	REPLACE 2
9 BALL TEMPERATURE OVERHEAT	TRANSISTOR	REPLACE 3

WIRE CONNECTION

A – POWER CORD

B – MAIN SWITCH


C – 9 BALL CABLE SOCKET

D – KEY PCB

NM300



KEY PCB



NOTE : ALWAYS CHECK THE VOLTS SWITCH BEFORE ASSEMBLY (MUST SET TO 220V)



PROBLEM	PARTS	SOLUTION
NO RESPONSE	TACK SWITCH	REPLACE BROKEN TACK SWITCH

NM300



9 BALL PROJECTOR

WIRE CONNECTION: SERIES TYPE →



9B CABLE CONNECTION TO PROJECTOR

1 – BROWN & BULB

2 – BLACK & SENSOR

3 – RED & SENSOR


4 – WHITE & BULB




9B CABLE CONNECTION TO PROJECTOR

PROBLEM	PARTS	SOLUTION
03 / 0E	BULB	REPLACE BULB
ON/OFF LIGHTS WHEN MOVING CABLE	9 BALL CABLE	REPAIR / REPLACE CABLE
2E	SENSOR	REPLACE SENSOR

Test if Bulb is working or not



Good – tester response
Bad – tester no response



Lamp lifespan in color



New



Good



Bad

THANK YOU!